


# XPLORE 1.1™

Physical Quantities Measurement System (PQMS)

User Manual



Precision-Quazar Tech Pvt. Ltd.  
New Delhi, India



# Contents

<b>1</b>	<b>Introduction to XPLORE 1.1<sup>TM</sup> System</b>	<b>1</b>
1.1	XPLORE 1.1 <sup>TM</sup> System Components . . . . .	2
1.2	Cryostat and Vacuum Manifold . . . . .	3
1.3	Inserts . . . . .	4
1.4	Power Supply Module (XPS) . . . . .	6
1.5	Temperature Controller Module (XTCON) . . . . .	7
1.6	Source & Measurement Module (XSMU) . . . . .	8
1.7	Lock-in Amplifier Module (XLIA) . . . . .	10
1.8	Electrical connections . . . . .	11
<b>2</b>	<b>4-Probe Resistance Measurement</b>	<b>15</b>
2.1	Mounting sample . . . . .	16
2.2	Acquiring $I$ - $V$ data . . . . .	17
2.3	Acquiring $R$ - $T$ data . . . . .	18
<b>3</b>	<b>Magnetic AC Susceptibility Measurement</b>	<b>21</b>
3.1	Mounting Sample . . . . .	22
3.2	Lock-in amplifier settings . . . . .	23
3.3	Acquiring $\chi$ - $T$ data . . . . .	24
<b>A</b>	<b>Installation requirements</b>	<b>27</b>
<b>B</b>	<b>Component checklist</b>	<b>29</b>
<b>C</b>	<b>System specification</b>	<b>31</b>
C.1	Cryostat and temperature controller . . . . .	32
C.2	$I$ - $V$ source and measurement unit for 4-probe $R$ - $T$ . . . . .	32
C.3	Lock-in amplifier for AC-susceptibility ( $\chi$ - $T$ ) measurement . . . . .	33
C.4	Data acquisition and control software: Qrius 1.1 . . . . .	34
C.5	Other user provided accessories required to run the system . . . . .	34
<b>D</b>	<b>Cable details</b>	<b>35</b>



# List of Figures

1.1	Major sub-units of XPLORE 1.1 <sup>TM</sup> setup . . . . .	2
1.2	Cryostat and vacuum manifold . . . . .	3
1.3	Four-probe I-V and R-T measurement insert . . . . .	4
1.4	AC susceptibility measurement insert . . . . .	5
1.5	XPLORE 1.1 <sup>TM</sup> power supply front & rear view . . . . .	6
1.6	XPLORE 1.1 <sup>TM</sup> temperature controller front & rear view . . . . .	7
1.7	XPLORE 1.1 <sup>TM</sup> source & measurement unit front & rear View . . . . .	8
1.8	XPLORE 1.1 <sup>TM</sup> lock-in amplifier front & rear view . . . . .	10
1.9	XPLORE 1.1 <sup>TM</sup> electronics power supply and computer interface connection	11
1.10	XPLORE 1.1 <sup>TM</sup> electronics to cryostat connection for R-T measurement .	12
1.11	XPLORE 1.1 <sup>TM</sup> electronics to cryostat connection for $\chi$ -T measurement .	13
2.1	Sample mount for $I$ - $V$ and $R$ - $T$ measurement . . . . .	16
2.2	Settings to acquire $I$ - $V$ data. . . . .	17
2.3	Settings to acquire $R$ - $T$ data. . . . .	18
3.1	Lock-in amplifier block diagram . . . . .	23
3.2	Launching Qrius $\chi$ -T module . . . . .	24



# Chapter 1

## Introduction to XPLORE 1.1<sup>TM</sup> System

## 1.1 XPLORE 1.1™ System Components




Figure 1.1: Major sub-units of XPLORE 1.1™ setup

Component	Description
<b>Mechanical</b>	
Cryostat	Non-magnetic double walled construction.
Vacuum manifold	To evacuate the cryostat to reduce evaporation of liquid N <sub>2</sub> .
<b>Electronics</b>	
Power supply	Powers up all XPLORE 1.1™ electronics modules.
Temperature controller	Controls temperature of cryostat heater and reads sample temperature.
I-V source & measurement unit	Measures resistance of the sample under constant voltage or constant current mode.
Lock-in amplifier	Measures output from secondary winding of magnetic AC susceptometer coil.



## 1.2 Cryostat and Vacuum Manifold




Figure 1.2: Cryostat and vacuum manifold

Item	Description	Port	Remark
A	Sample insertion port.	KF-25	Supplied
B	Heater chamber evacuation port.	KF-16	Supplied
C	Sample chamber evacuation port.	KF-16	Not supplied
D	Heater electrical connection feed-through.	8-pin BNC	Supplied
E	Cryostat stand.		
F	Non-magnetic double-walled cryostat.		
G	Pirani gauge mounting port.	KF-10	Supplied
H	Release valve.	KF-10	Supplied
I	Vacuum pump connection port.	KF-25	Not supplied
J	Butterfly valve.		
K	Cryostat vacuum hose connection port.	KF-25	Supplied

### 1.3 Inserts




Figure 1.3: Four-probe I-V and R-T measurement insert




Figure 1.4: AC susceptibility measurement insert

## 1.4 Power Supply Module (XPS)




Figure 1.5: XPLORE 1.1™ power supply front & rear view

## 1.5 Temperature Controller Module (XTCON)




Figure 1.6: XPLORE 1.1™ temperature controller front & rear view

Parameter	Value
<b>Temperature measurement and control</b>	
Temperature range	77 K – 450 K
Sensors	Pt100 and K-type thermocouple
Temperature resolution	0.01 K
Temperature stability	Better than $\pm 0.1$ K in isothermal mode
Temperature ramp-rate	1 K/min – 10 K/min in steps of 0.1 K/min
Control algorithm	PID
PID parameters	User configurable
<b>Heater</b>	
Heater power	40 Watt
<b>Others</b>	
Computer connectivity	USB

## 1.6 Source & Measurement Module (XSMU)




Figure 1.7: XPLORE 1.1™ source & measurement unit front & rear View

Parameter	Value
Source mode	Current or voltage
<b>Current source specification</b>	
Current source range(s)	100 $\mu$ A, 1 mA, 10 mA
Current set-point resolution	Better than 0.05 % of full-scale
Voltage compliance	10V
<b>Voltage source specification</b>	
Voltage source range(s)	10 V
Voltage set-point resolution	Better than 0.005 % of full-scale
Current compliance	10 mA
<b>Ammeter specification</b>	
Current measurement range(s)	100 $\mu$ A, 1 mA, 10 mA
Current measurement resolution	6½ digit
<b>Voltmeter specification</b>	
Voltage measurement range(s)	100 mV, 1 V, 10 V

I-V source and measurement unit specification continued on next page ...

---

... I-V source and measurement unit specification continued from previous page

---

Parameter	Value
Voltage measurement resolution	6½ digit
Input impedance	$> 10^{12} \Omega$
<b>Ohmmeter specification</b>	
Resistance measurement range(s)	10 $\Omega$ , 100 $\Omega$ , 1 k $\Omega$ , 10 k $\Omega$ , 100 k $\Omega$ , 1 M $\Omega$ , 10 M $\Omega$
Resistance measurement resolution	6½ digit
<b>Others</b>	
Computer connectivity	USB

---

## 1.7 Lock-in Amplifier Module (XLIA)




Figure 1.8: XPLORE 1.1™ lock-in amplifier front & rear view

Parameter	Value
<b>Primary driver</b>	
Frequency range(s)	10 Hz – 10 KHz
Reference output	2V peak-to-peak
Reference generation	Direct digital synthesis
<b>Lock-in detection</b>	
Integration time	1 sec
Pre-amp gain(s)	1, 10, 100
Post-amp gain(s)	1, 10, 100
<b>Measurement</b>	
Measurement resolution	5½ digit
Measured quantities	In phase, quadrature phase, amplitude, and phase.
<b>Others</b>	
Computer connectivity	USB



## 1.8 Electrical connections




Figure 1.9: XPLORE 1.1<sup>TM</sup> electronics power supply and computer interface connection

Figure no.	Description
Figure 1.9	Shows basic power supply and computer interface connection irrespective of the nature of the experiment.
Figure 1.10	Shows connections between electronics modules, and cryostat and inserts for $I$ - $V$ and $R$ - $T$ experiments.
Figure 1.11	Shows connections between electronics modules, and cryostat and inserts for $\chi$ - $T$ experiment.




Figure 1.10: XPLORE 1.1™ electronics to cryostat connection for R-T measurement




Figure 1.11: XPLORE 1.1™ electronics to cryostat connection for  $\chi$ -T measurement



# Chapter 2

## 4-Probe Resistance Measurement

## 2.1 Mounting sample




Figure 2.1: Sample mount for  $I$ - $V$  and  $R$ - $T$  measurement

- Paste the sample on the sample mount as shown in Figure 2.1. Position it as close to the thermocouple as possible to avoid lag in temperature readout.
- Solder four thin wires to the contact pads  $I+$ ,  $I-$ ,  $V+$ , and  $V-$ .
- Connect their other ends to the sample using silver epoxy, indium solder or any other technique, which ensures a good electrical contact of the wires to the sample.
- Gently slide down the four-probe measurement insert in the cryostat and clamp it using its top flange (KF-25).
- Connect the insert with the XSMU and XTCON as described in Section 1.8.

## 2.2 Acquiring $I$ - $V$ data




Figure 2.2: Settings to acquire  $I$ - $V$  data.

- Launch Qrius and open **4-probe  $I$ - $V$  and  $R$ - $T$**  module as shown in Figure 2.2. Consider the dotted rectangle.
- In the **Current Ramp Settings** frame on the four-probe window, enter the initial and final values of the current with the step increment in current value.
- Select the **V Measurement Full Scale Range** in order to ensure that the measurement does not saturate during the sweep.
- Start  $I$ - $V$  data logging by pressing the **Start Current Ramp** button.
- Check if the data obtained shows the expected curve.
- If the curve is a flat-line (saturated measurement inputs), then increment the **V Measurement Full** setting.
- Restart the  $I$ - $V$  log by again pressing the **Start Current Ramp** button.

## 2.3 Acquiring $R$ - $T$ data

**Note:** Before taking an  $R$ - $T$  plot, log  $I$ - $V$  data to ensure that the contacts are fine, and current and voltmeter settings are within range.




Figure 2.3: Settings to acquire  $R$ - $T$  data.

- Launch Qrius and open **4-probe  $I$ - $V$  and  $R$ - $T$**  module as shown in Figure 2.3. Consider the solid rectangles.
- Enter the current to be supplied to the sample in the **Set Current** entry box & Press **return** key. The box turns yellow to indicate that the specified current has been set.
- Select the (optional) **V Measurement Full Scale Range**.
- Enter the **final temperature** and **ramp rate** values on the XTCON interface. Press **return** key to transmit these settings to the XTCON hardware.
- On the **four-probe interface**, tick on the box which directs it to use the temperature ramp settings provided in the XTCON interface.
- Click on the **Start  $R$ - $T$  Plot** button to initiate  $R$ - $T$  data logging.



- The final temperature and ramp rate can be modified on-the-fly during an experiment.
- To stop  $R$ - $T$  data logging, click on Stop  $R$ - $T$  Plot button.

**Please note**

- A negative ramp can also be entered as long as the final temperature is lower than the current temperature.
- The initial temperature setting is ignored. Current temperature is used as the initial temperature in a ramp.
- Apply a lower ramp rate ( $< 2^{\circ}\text{C}/\text{min}$ ) to avoid temperature sensor lag.



# Chapter 3

## Magnetic AC Susceptibility Measurement

### 3.1 Mounting Sample

- Encapsulate the sample in a butter-paper wrapping and insert it into the sample holder pipe sleeve.
- Ensure that the sample position in the sleeve is such that it can traverse freely on either side, from the middle of the primary coil.
- See to it that no other magnetic material is introduced in the coil except the test sample.
- Slowly lower the sample into the coil while observing the amplitude output of the XLIA module on its LCD display.
- Fix the position of the sample when the maximum value on the amplitude is observed.
- Position the cryostat and the susceptibility insert in a position that the sample location remains undisturbed throughout the experiment.

## 3.2 Lock-in amplifier settings




Figure 3.1: Lock-in amplifier block diagram

### 3.3 Acquiring $\chi$ -T data




Figure 3.2: Launching Qrius  $\chi$ -T module

- Launch Qrius and open **Magnetic AC susceptibility ( $\chi$ -T)** module as shown in Figure 3.2. Consider the solid rectangles.
- Enter the **final temperature** and **ramp rate** values on the XTCON interface. Press **return** key to transmit these settings to the XTCON hardware.
- Set lock-in amplifier parameters, as shown in Figure 3.1.
- On the **AC susceptance** window, tick on the box which directs it to use the temperature ramp settings provided in the XTCON interface.
- Click on the **Start  $\chi$ -T Plot** button to initiate data logging.
- Click on the **Stop  $\chi$ -T Plot** to terminate an ongoing run.





# Appendix A

## Installation requirements

---

### Power supply

- 1 8 nos. 240 VAC 5 A 3-pin sockets.
- 2 Ground should not be more than 5 V w.r.t neutral.
- 3 3 KVA UPS back-up of all sockets if power failure is frequent.

### Space

- 1 One 5 ft × 3 ft dedicated table for comfortable accommodation of electronics, vacuum manifold, and computer.

### Computer Requirements

- 1 DVD drive
- 2 2 GB RAM
- 3 3 free USB ports
- 4 30 GB free partition on hard disk
- 5 1024x768 (or better) monitor
- 6 We would recommend a dedicated PC for XPLORE 1.1 system.

### Additional Requirements

- 1 3 liter Liquid N<sub>2</sub> for experiments where cooling is required.
  - 2 Soldering iron, GE-varnish, silver paste, laminated copper wire, etc. for mounting samples.
-



Appendix

B

## Component checklist

### Electronics boxes

1	Power supply	1
2	Temperature controller	1
3	I-V source & measurement unit	1
4	Lock-in amplifier	1
5	Pirani gauge readout unit	1

### Electrical cables

1	Power cord	4
2	Cryostat heater cable	1
3	4-probe I-V and R-T insert cable	1
4	AC susceptometer coil insert cable	1
5	AC susceptometer sample insert cable	1
6	Power supply inter-connects	2
7	Pirani gauge inter-connect	1
8	USB cable	3

### Cryostat, inserts, and vacuum accessories

1	Cryostat	1
2	Cryostat stand	1
3	Vacuum manifold	1
4	Vacuum manifold base	1
5	Cryostat vacuum hose	1
6	Pirani gauge	1
7	Mounting base	1
8	C-clamp (3 inch)	2

---

... XPLORE 1.1™ component checklist continued from previous page

---

9	KF-10 clamp and center-ring	2
10	KF-16 clamp and center-ring	4
11	KF-25 clamp and center-ring	2
12	4-probe I-V and R-T insert	1
13	AC susceptometer coil insert	1
14	AC susceptometer sample insert	1
15	Liquid nitrogen dewar	1

#### Sample toolkit

1	Tool box	1
2	Kapton tape reel	1
3	Teflon tape reel	1
4	size-4 brush	1
5	GE-varnish	5 ml
6	Silver paste	2 ml
7	Copper wire	5 m
8	Screw driver set	1
9	Scissor	1
10	Wire cutter	1
11	Tweezers	1
12	Test samples	5
13	AC susceptometer sample holder sleeve	10

#### Vacuum pump and accessories

1	Vacuum pump	1
2	Bellow 1 meter, KF-25	1
3	KF-25 clamp and center-ring	2

---

# Appendix C

## System specification

## C.1 Cryostat and temperature controller

Parameter	Value
<b>Sample chamber</b>	
Dimension	21 mm (dia) × 100 mm (length), top loading
Construction	Non-magnetic, double walled
<b>Temperature range</b>	
Temperature range	80 K – 450 K
Cryogen	Liquid N <sub>2</sub>
Exchange gas	Air, port available for evacuation/flushing
<b>Temperature measurement and control</b>	
Sensors	Pt100 and K-type thermocouple
Temperature resolution	0.01 K
Temperature stability	Better than ±0.1 K in isothermal mode
Temperature ramp-rate	1 K/min – 10 K/min in steps of 0.1 K/min
Control algorithm	PID
PID parameters	User configurable
<b>Heater</b>	
Heater power	40 Watt
Winding type	Non-magnetic
<b>Others</b>	
Computer connectivity	USB

## C.2 I-V source and measurement unit for 4-probe $R$ - $T$

Parameter	Value
Source mode	Current or voltage
<b>Current source specification</b>	
Current source range(s)	100 $\mu$ A, 1 mA, 10 mA
Current set-point resolution	Better than 0.05 % of full-scale
Voltage compliance	10V

I-V source and measurement unit specification continued on next page ...

---

... I-V source and measurement unit specification continued from previous page

---

Parameter	Value
<b>Voltage source specification</b>	
Voltage source range(s)	10 V
Voltage set-point resolution	Better than 0.005 % of full-scale
Current compliance	10 mA
<b>Ammeter specification</b>	
Current measurement range(s)	100 $\mu$ A, 1 mA, 10 mA
Current measurement resolution	6½ digit
<b>Voltmeter specification</b>	
Voltage measurement range(s)	100 mV, 1 V, 10 V
Voltage measurement resolution	6½ digit
Input impedance	$> 10^{12} \Omega$
<b>Ohmmeter specification</b>	
Resistance measurement range(s)	10 $\Omega$ , 100 $\Omega$ , 1 k $\Omega$ , 10 k $\Omega$ , 100 k $\Omega$ , 1 M $\Omega$ , 10 M $\Omega$
Resistance measurement resolution	6½ digit
<b>Others</b>	
Computer connectivity	USB

---

### C.3 Lock-in amplifier for AC-susceptibility ( $\chi$ - $T$ ) measurement

Parameter	Value
<b>Sample</b>	
Sample size	5 mm (dia) $\times$ 10 mm (length)
Field strength	4 Oe RMS field at 80 Hz
<b>Primary driver</b>	
Frequency range(s)	10 Hz – 10 KHz
Reference output	2V peak-to-peak
Reference generation	Direct digital synthesis

---

AC-susceptometer specification continued on next page ...

---

---

... AC-susceptometer specification continued from previous page

---

Parameter	Value
<b>Lock-in detection</b>	
Integration time	1 sec
Pre-amp gain(s)	1, 10, 100
Post-amp gain(s)	1, 10, 100
<b>Measurement</b>	
Measurement resolution	16-bit
Measured quantities	Both amplitude and phase
<b>Others</b>	
Computer connectivity	USB

---

## C.4 Data acquisition and control software: Qrius 1.1

Computer automation of  $R$ - $T$  and  $\chi$ - $T$  measurements using above mentioned modules have been provided by Qrius 1.1 software. It allows easy control of all experimental parameters and real time recording and plot of physical quantities. The software runs on Ubuntu 12.04.

## C.5 Other user provided accessories required to run the system

The following is a list of accessories that *need to be provided by the user* for proper operation of XPLORE-1.1 PQMS setup.

- Rotary vacuum pump.
- Minimum of 3 liters of liquid nitrogen.
- Computer with DVD drive, USB 2.0 port, 2GB RAM, 100GB of hard-disk space and Pentium dual core or above processor.



# Appendix D

## Cable details

